



EZCOOL

INSTRUCTIONS OF USE

BEFORE PRF PREPARATION

- 1) Fill the EZCOOL cup with 2/3 cup of cold water.
- 2) 20 min prior to the treatment, turn on the EZCOOL device by pressing the ON/OFF touch sensor on the front. When on, a blue light will appear.



- 3) Once the device is on, proceed with the blood draw and centrifugation as usual.

DURING PRF PREPARATION

- 1) While preparing the PRF, the EZCOOL device will chill the water to the precise preset temperature.

AFTER PRF PREPARATION

- 1) Once centrifugation time is up (+/- 5 min), place your EZPRF tubes in the pre-chilled EZCOOL water cup. The PRF will stay in its liquid form for 40-60 minutes after centrifugation.
- 2) Once ready to use, take the EZPRF tube out of the EZCOOL and withdraw the PRF.

A PRODUCT OF



COSMOFRANCE